


# GENERIC PROTOCOL :

## Electrical Stimulation for Chronic Wounds

### MONOPOLAR ELECTRODE SET UP



**Active:** Smaller electrode placed directly in/over ulcer  
composed of sterile/single use conductive material

Types:

- saline soaked gauze, hydrogel gauze, hydrofibre
- covered with self adhesive or carbon electrode
- silver nylon dressing

**Dispersive:**

LARGE that is at least twice the size of active  
placed on intact healthy skin  
proximal to wound at least 20 cm away

Types:

- self adhesive (5" x 8")
- carbon coated electrode over damp cloth

